



StabilZyme® AP Conjugate Stabilizer-CMIT/MIT

Product insert

Product number:

SA01-CF01-0050 (50 mL)
SA01-CF01-1000 (1000 mL)
SA01-CF01-2000 (2000 mL)

Intended use:

StabilZyme® conjugate/protein stabilizers protect the entire conjugate by preventing the loss of catalytic activity and maintaining the structural integrity of the protein in solution. These stabilizers allow for the storage of conjugated proteins at lower use concentrations, extended shelf life and increased signal-to-noise ratios for improved assay performance.

The StabilZyme AP Conjugate Stabilizer-CMIT/MIT is optimized to maintain the conformation of antibodies/antigens conjugated to alkaline phosphatase (AP) in solution. Adding the conjugate directly into the AP Conjugate Stabilizer allows for stable storage of the conjugate at a working concentration and eliminates the need for subsequent dilutions.

Product stability, storage and specifications:

Product stability	Stable for 2.5 years from date of manufacture
Storage	Product should be stored at 2-8°C or room temperature
Specification	Bovine Protein: Contains bovine serum albumin Product Buffer: MES pH: 5.5 - 6.5 Preservative: 0.00144% CMIT/MIT
Notes	Please note that SurModics Stabilization Products are shipped to customers at ambient temperature. Extensive stability studies have shown that prolonged storage at ambient temperature will not affect the product quality or efficacy.

Recommendations for use:

Aseptically pour off the desired volumes of material needed for the application and allow the product to equilibrate at room temperature prior to use.

The following are general guidelines only.

- 1) Use StabilZyme AP Conjugate stabilizer-CMIT/MIT at 100% concentration for optimal stability and performance.
 - If necessary, adjust the pH of the solution to align with the properties of the protein for optimal performance.
- 2) Dilute the AP conjugate/protein/antibody coated particles to a working concentration in the StabilZyme AP Conjugate Stabilizer-CMIT/MIT solution.
- 3) Use the diluted conjugate solution according to the lab defined protocol for the assay.
- 4) For optimal performance, store the conjugate solution at 2-8°C and protect from direct exposure to light.

Additional considerations:

An accelerated stability test is recommended to adequately assess the stability of the conjugate in solution. Samples are stored at 37°C and 4°C for up to 90 days (longer if you prefer). At set intervals, the percent



StabilZyme® AP Conjugate
Stabilizer-CMIT/MIT
Product insert

retained activity is calculated by dividing the 37°C sample optical density (O.D.) by the 4°C sample O.D. and plotted over time.

Related products:

In-Solution Protein Stabilizers & Diluents:
StabilZyme® Protein Free AP Stabilizer (SAPF)
Assay Diluent-HAMA Blocker-CMIT/MIT (SM01-CF02)
Blockers/Stabilizers:
StabilGuard® Immunoassay Stabilizer BSA-Free (SG01)
StabilCoat® Immunoassay Stabilizer (SC01)
Substrates:
BioFX® pNPP Yellow One Component Microwell Substrate with Stabilizing Pellets (PNPS)
BioFX® Stabilizing Pellets for PNPS (STAB)
BioFX® BCIP/NBT Blue One Component AP Membrane Substrate (BCIB)
BioFX® BCIP/NBT Purple One Component AP Membrane Substrate (BCID)
BioFX® Chemiluminescent Supersensitive AP Microwell/Membrane Substrate-450 nm (APS4)
BioFX® Chemiluminescent Supersensitive AP Microwell/Membrane Substrate-540 nm (APS5)
BioFX® Chemiluminescent Ultra-Sensitive AP Microwell/Membrane Substrate-450 nm (APU4)
BioFX® Chemiluminescent Ultra-Sensitive AP Microwell/Membrane Substrate-450 nm (APU5)

Products distributed by



DIARECT AG
Bötzing Str. 29 B
79111 Freiburg
Germany

Tel. +49 761 47979-0
Fax +49 761 47979-29
www.diarect.com
orders@diarect.com
info@diarect.com

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made according to DIARECT's Terms and Conditions of Sale.